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User Manual / Safety Data Sheet
www.microceram.de

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GB:
 The MicroCeram Dentalblanks are Yttria stabilized, pre-sintered Zirconium oxide, suitable for manufacturing wide variety of dental restoratives viz. copings, crowns, bridges and other frameworks. The Company are certified, by TÜV Rheinland LGA Products GmbH and the Blanks are compliance with following ISO Standards: ISO 13356 Second Edition 2008-06-01 Implants for surgery-ceramics material based on Yttrium stabilized tetragonal Zirconium - ISO 6872:2009-01 Dentistry-Ceramic materials.

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GB: Handling

MicroCeram dental blanks can be milled on any manually operated copy milling systems as well as on dental CAD/CAM, to produce single copings, crowns, bridges and other frameworks.

It is imperative that the blanks be checked to ensure that their optical appearance is flawless. If a blank which has been damaged (during transport) is used, a subsequent complaint will no longer be accepted.

The shrinkage factor of MicroCeram blanks is given to the batch. This has to be regulated on the milling system or input into the CAD software of the CAD/CAM system.

According to requirements, the frameworks can be dyed with Liquid dyeing solutions, either partially or completely (Cf. processing instructions of the dyeing solutions).

The advantages of MicroCeram blanks are:

- Rawmaterial quality
- Handling
- High bendingstrength data of Zirconium
- Good biological compatibility of the product
- Resistance to corrosion
- Its light color and the individual coloring of the blocks
- Its translucency with thin wall thicknesses

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Contraindication

- Insufficient oral hygiene
- Insufficient preparation results
- Insufficient tooth structure
- Insufficient space available
- Bruxism
- Known intolerance to the constituents
- Care with use during pregnancy or in children

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GB: Preparation / Material-specific framework parameters

- The preparation must be performed with either a chamfer or a shoulder with rounded internal angle.
- The vertical preparation angle should be at least 3°.
- All transitions from the axial to the occlusal or incisal areas must be rounded off. Flat or plane surfaces are advantageous.

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GB: Sinterprocess

The sintering of the ceramics frameworks take place according to the following heating program: From room temperature 20 °C to final sintering temperature 1450 °C - Heating rate of the kiln 4-8 K/min final temperature - holding time 2 hours - Kiln cooling down time to temperature < 200 °C, approx 5 °C / min

Lay the frameworks on loose high purity (99.7 % Al₂O₃) ceramic micro beads loaded in refractory ceramic sagger. Apply gentle pressure on the frameworks so that they take proper support on the loose ceramic beads. Ensure that no ceramic beads get jammed in the internal frameworks cavities. (Please read the instructions of the sintering)

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GB: Post processing / Fine working

After milling with diamond tools, the frames may be manually ground if necessary to add more finish to the shape. The framework must always be wet ground to avoid overheating. After checking fitting of frameworks accuracy, contact points an occlusion, clean the interior cavity of ceramic framework by the corundium or suitable jets (Aluminium oxide 110 µm / pressure ≤ 2 bar) and subsequent operation / cleaning of the same.

Porcelain bonding / Veneer

All known standard veneer ceramics for zirconium oxide can be used as long as their thermal expansion coefficient is between 10 and 10.5. Clean the framework after final shaping carefully with hot steam. Cover the zirconium frame completely with a layer of ceramic.

GB: Fastening instructions

The Restoration can be fastened non-adhesively with glasionomer or zincphosphate cements, or adhesively with the self-curing composite. In all this case the manufacturer's processing instructions must be observed without fail.

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GB: Warning and safety instructions

The dental blanks must be used according to instructions. We accept no liability for damages resulting from improper handling or processing. The user is obliged to check the product before use. Furthermore, our liability for the accuracy of this information is independent of reason and, if legally permissible, is limited in each case the value of the invoiced without VAT. We are not liable for lost profits, indirect damages, consequential damages or claims by third parties against the buyer. Damages based on claims are given only in cases of intent or gross negligence. When processing, please use personal protective clothing (dust mask, protective goggles).

Storage

Store in a clean and dry place (+5 °C and 50 °C) in its original packaging.

Shelf life

Shelf life of up to 5 years from the date of manufacture.

The current "User Manual" you can be downloaded at: www.microceram.de

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We reserve the right to make any alterations which may be required due to technical improvements.